

CADD Prizm[®] VIP Ambulatory Infusion Pump

Models 6100 and 6101

Intermittent Delivery Mode

1. Reservoir Volume
2. Dose Volume
3. Dose Duration
4. Dose Cycle
5. KVO Rate
6. Next Dose Start Time
7. Dose Rate (ml/hr)
8. Milliliters Given

TPN Delivery Mode

1. Reservoir Volume
2. Infusion Volume
3. Infusion Period
4. Taper-Up Period
5. Taper-Down Period
6. Plateau Rate
KVO Rate
7. Milliliters Given

Continuous Delivery Mode

1. Reservoir Volume
2. Continuous Rate
3. Milliliters Given

PCA Delivery Mode

1. Reservoir Volume
2. Units
3. Concentration
4. Continuous Rate
5. Demand Dose
6. Demand Dose Lockout
7. Demand Doses Per Hour
8. Dose Counters
Given/Attempted
9. (Units) Given
10. New Patient Marker
(If Extended History is ON)

Deltec

Deltec, Inc.
1265 Grey Fox Rd., St. Paul, MN 55112 USA

1-800-426-2448 www.deltec.com
C0200

CADD-Prizm design logo and CADD-Prizm are Deltec, Inc. trademarks.
©2001 Deltec, Inc. All rights reserved. Printed in U.S.A. Revised 12/01

smiths
Smiths Medical – a part of Smiths Group plc